



PRODUCT CATALOG



ARC Medical, Inc.

circuitGuard™ ♦ **ThermoFlo™**

4296 Cowan Road Tucker, GA 30079

Order Number: 800.950.ARC1 (2721)

General Information Voice Mail: 800.950.2720 Ext 727

E-Mail: arcinfo@arcmedical.com

Web: www.arcmedical.com



PRODUCTS FOR RESPIRATORY CARE

THERMOFLO™ FILTER



Maintains sputum consistency and provides therapeutic humidification in a filtered device. Zero ET tube occlusions reported.

ORDER INFORMATION:

#6000	THERMOFLO FILTER W/HCH	40/CASE
#6001	THERMOFLO FILTER W/HCH EcoPk	100/CASE
#6001V	THERMOFLO FILTER W/HCH EcoPk	200/CASE
#6020	THERMOFLO FILTER W/ETCO ₂ GS PORT	40/CASE

Specifications at high V_E:

Moisture Returned H₂O/L: 33 mg @ 10L V_E & 31 mg @ 20L V_E

Resistance: @ 30 LPM - 1.2 cm

Dead Space: 75 ml Weight: 33 grams Housing: Clear

Tidal Volume Range: 250 ml - 1500 ml

Bacterial/Viral Efficiency: 99.9%+ TB Filtration Efficiency: 99.9%+



THERMOFLO™ FILTER SPECIAL

The "S" Version provides a collection chamber to help guard against increasing resistance due to secretions. Zero ET tube occlusions reported.

ORDER INFORMATION:

#6000S	THERMOFLO SPECIAL	40/CASES
#6001S	THERMOFLO SPECIAL EcoPAK	100/CASES
#6000SL	THERMOFLO SPECIAL W/LUER PORT	40/CASES

Specifications at high V_E:

Moisture Returned H₂O/L: 33mg @ 10L V_E & 31mg @ 20L V_E

Resistance: @ 30 LPM - 1.2 cm

Dead Space: 75 ml Weight: 33 grams Housing: Clear

Tidal Volume Range: 250ml - 1500ml

Bacterial/Viral Efficiency: 99.9%+ TB Filtration Efficiency: 99.9%+



THERMOFLO™ 1 FILTER & HCH

The small size and low profile are excellent for pediatric, adult and, anesthesia patients.

ORDER INFORMATION:

#6100	THERMOFLO1 FILTER HCH	50/CASE
#6101	THERMOFLO1 FILTER HCH EcoPk	100/CASE
#6120	THERMOFLO1 FILTER HCH ETCO ₂	50/CASE
#6121	THERMOFLO1 FILTER HCH ETCO ₂	100/CASE

Specifications at high V_E:

Moisture Returned H₂O/L: 30.2 mg @ 10L V_E

Resistance: @ 30 LPM - 1.2 cm

Dead Space: 30 ml Weight: 14 grams Housing: Clear

Tidal Volume Range: 100 ml - 1200 ml

Bacterial/Viral Efficiency: 99.9%+ TB Filtration Efficiency: 99.9%+



PRODUCTS FOR RESPIRATORY CARE



THERMOFLO™ HCH (NO FILTER)

Therapeutic humidification in an economically priced HCH. Useful when a filter is not needed.

ORDER INFORMATION:

#6070	THERMOFLO HCH	40/CASE
#6071	THERMOFLO HCH EcoPAK	100/CASE
#6072	THERMOFLO HCH W/ETCO ₂ PORT	40/CASE
#6073	THERMOFLO HCH W/ETCO ₂ PORT	100/CASE

Specifications at high V_E:

Moisture Returned H₂O/L: 30.2mg @ 10L V_E

Resistance: @ 30 LPM - 0.4 cm

Dead Space: 75 ml

Weight: 32 grams

Housing: Clear

Tidal Volume Range: 250 ml - 1500 ml



THERMOFLO™ 1 HCH (NO FILTER)

The small size and low profile are excellent for pediatric and adult patients. Economically priced and useful when a filter is not needed.

ORDER INFORMATION:

#6060	THERMOFLO 1	50/CASE
#6061	THERMOFLO 1	100/CASE
#6160	THERMOFLO 1 W/ ETCO ₂ PORT	50/CASE
#6161	THERMOFLO 1 W/ETCO ₂ PORT	100/CASE

Specifications at high V_E:

Moisture Returned H₂O/L: 33mg @ 10L V_E & 31mg @ 20L V_E

Resistance: @ 30 LPM - 0.7cm

Dead Space: 30ml

Weight: 14 grams

Housing: Clear

Tidal Volume Range: 100 ml - 1200 ml



THERMOFLO™ NEO HCH

Transport moisture for patients who need it the most. Unexcelled low resistance to flow, critical in the care of special care babies. Non-absorbing HCH material reduces resistance problems. .

ORDER INFORMATION:

#6215 THERMOFLO NEO

#6220 THERMOFLO NEO W/ETCO₂ PORT

BOTH ARE PACKAGED 30/CASE

Specifications:

Moisture @ V_T 20ml = 32.5mg H₂O/L air @/V_T 25ml = 31.0 H₂O/L air


Resistance @ 7.5 L/min = 0.35cm H₂O @ 15L/min 0.95cm H₂O

Dead Space: 2.5 ml

Weight: 3 grams

Housing: Clear

Tidal Volume Range: 5ml - 50 ml



THERMOFLO™ TRACH HCH

Available in either a **basic trach** or **CRV** (cough relief valve) trach. The CRV may be used for suctioning without removing the device, and features an optional integrated O₂ adapter for supplemental O₂ delivery. The HCH material is non-absorbing, helping to eliminate resistance. These are designed for spontaneously breathing patients through a tracheotomy tube.

THERMOFLO™ TRACH BASIC:

#6240 PACKED 30/BOX #6240E 240 UNITS (PACKED 8 BAGS OF 30 EACH)

THERMOFLO™ TRACH BASIC OXYGEN ADAPTER:

#6258 25/BOX

THERMOFLO™ TRACH CRV:

#6241 WITH O₂ PORT 30/BOX; #6241E WITH O₂ PORT 240/BOX

#6242 WITHOUT O₂ PORT 30/BOX; #6242E WITHOUT O₂ PORT 240/BOX

Specifications:

Humidity V_T 500ml: 26mg H₂O / L air
Resistance @ 0.7 cm H₂O @30 LPM
Dead Space: 12 ml

TAPER DRAPE FLEX TUBE

ORDER #5010T 50/CASE





PRODUCTS FOR ANESTHESIA



CIRCUITGUARD™ HCH & FILTER

Protects patients and anesthesia equipment from contamination. The three main components of circuitGuard™ are manufactured as one unit using our exclusive sure-fit system. No glue or welding is needed, eliminating possible fail points found on other products. This particular model has a HCH.

ORDER INFORMATION:

#6126 CIRCUITGUARD STRAIGHT FILTER, HCH, ETCO₂ PORT & ELBOW 100/CASE

#6131 CIRCUITGUARD 90° RIGHT ANGLE FILTER, HCH & ETCO₂ PORT 100/CASE



Specifications:

Humidity Returned mg H₂O/L: 30.2mg @ 10L V_E

Resistance: @ 30 LPM - 1.2cm

Dead Space: 30ml

Weight: 14 grams

Housing: Clear

Tidal Volume Range: 100ml - 1200ml

Bacterial/Viral Efficiency: 99.9+% TB Efficiency 99.9+%



CIRCUITGUARD™ FILTER ONLY

Protects patients and anesthesia equipment from contamination. The three main components of circuitGuard™ are manufactured as one unit using our exclusive sure-fit system. No glue or welding is needed, eliminating possible failure points.

ORDER INFORMATION:

#7056 CIRCUITGUARD STRAIGHT FILTER ONLY, ETCO₂ PORT & ELBOW 100/CASE

Specifications:

Resistance: @ 30 LPM - 0.8 cm

Dead Space: 30ml

Weight: 13 grams

Housing: Clear

Tidal Volume Range: 100ml - 1200ml

Bacterial/Viral Efficiency: 99.9+% TB Efficiency 99.9+%



ANESTHESIA CIRCUITS

Our anesthesia circuits are available in your choice of expandable or corrugated tubes. These circuits have no ports, no elbow and no filter. They are packaged with the nrye piece, attached tubes with cuff adapters and latex free re-breathing bag.

Expandable Circuits		
X40LF3	30 case	40"flex, 3L latex free
X72LF3	30 case	72"flex, 3L latex free
X96LF3	30 case	96"flex, 3L latex free
Corrugated Circuits		
D40LF3	30 case	40"hose, 3L latex free
D60LF3	30 case	60"hose, 3L latex free
D72LF3	30 case	72"hose, 3L latex free

SOFT AIR CUSHION FACE MASKS WITH INFLATION VALVE



May be packaged with circuitGuard™ for your convenience.

ORDER INFORMATION: PACKAGED 50/CASE

- #5001 = SIZE 1**
- #5002 = SIZE 2**
- #5003 = SIZE 3**
- #5004 = SIZE 4**
- #5005 = SIZE 5**

NOTES

To Order Call:
800.950.ARC1 (2721)



ARC Medical, Inc.

4296 Cowan Road
Tucker, GA 30084

General Info Voice Mail: 800.950.2720 Ext 272

Order: 800.950.ARC1 (2721)

E-Mail: arcinfo@arcmedical.com

www.arcmedical.com