




# ThermoFlo™ TrachCRV

Now with Cough Relief Valve

**Always Low  
Prices!**

ICU / ENT / Home Care

**ThermoFlo™ Trach CRV**  
Actual Size



**ThermoFlo™ Trach CRV O<sub>2</sub>**  
Integral Adapter for Supplemental O<sub>2</sub>

- O<sub>2</sub> Humidified before delivery to patient
- Inspiratory Gases are Optimally Humidified
- High Moisture Output with Low Resistance to Flow
- Adult and Pediatric Patients
- Dead Space: 12 ml
- Humidification V<sub>T</sub> 500ml: 26mg H<sub>2</sub>O / L air
- Resistance: 0.7cm H<sub>2</sub>O @ 30 LPM

<u>Description</u>	<u>Box of 30</u>	<u>8 boxes of 30</u>
ThermoFlo™ Trach Basic	6240	6240E
ThermoFlo™ Trach CRV with O <sub>2</sub>	6241	6241E
ThermoFlo™ Trach CRV	6242	6242E

*Humidification @ Lowest Daily Cost*

Orders: 800-950 2721  
Phone: 404-373 8311  
Fax: 404-373 8385  
Web: [arcmedical.com](http://arcmedical.com)

**ARC Medical, Inc.**  
4296 Cowan Road  
Tucker, GA 30084